
SPIRALES *FREE PLATINUM SPIRALES*





SPIRALES

BALT's spirales are composed of a fine platinum wire. They are designed for intravascular use in neuroradiological of peripheral zones, for the treatment of arterio-vascular malformations and aneurysms. Three ranges of coils are available: free platinum spirales with or without fibers, and Flow Coils.

SPI & CIRRUS Ranges

The free radiopaque spirales, with or without fibers, are put on a straightening mandrel enabling an easy introduction into a catheter.

SPI Range (without fibers)

SPI.037" -	5F	catheters	(1,1-1	,25mm)	– GW	.038"

Reference	PrimarySpring Ø	Coil Ø mm	Coil length cm
SPI 5X300P37		5	30
SPI 6X30P37		6	3
SPI 6X60P37		6	6
SPI 6X120P37		6	12
SPI 6X300P37		6	30
SPI 8X60P37	0.82 mm / [.032"]	8	6
SPI 8X120P37	0.02 11117 [.002]	8	12
SPI 8X300P37		8	30
SPI 10X60P37		10	6
SPI 10X120P37		10	12
SPI 10X300P37		10	30
SPI 10X400P37		10	40
SPI 16X60P37		16	6
SPI 16X300P37		16	30
SPI 16X450P37		16	45

SPI.018" - VASCO+18 - GW .014"

Reference	PrimarySpring Ø	Coil Ø mm	Coil length cm
SPI 2X25P18	0.36 mm / [.014"]	2	2,5
SPI 2X50P18		2	5
SPI 3X50P18		3	5
SPI 3X80P18		3	8
SPI 4X120P18		4	12
SPI 5X150P18	0.40 mm / [.016"]	5	15
SPI 7X150P18		7	15
SPI 9X150P18		9	15
SPI 11X150P18		11	15



SPI (without fibers)

CIRRUS Range (with fibers)

FIB.015" – 3F catheters (0,6-0,7mm) – GW .021"				
Reference	PrimarySpring Ø	Coil Ø	Coil length cm	
FIB 2X25P15	0.37 mm / [.015"]	2	2,5	
FIB 3X25P15		3	2,5	
FIB 3X50P15		3	5	
FIB 4X30P15		4	3	
FIB 4X50P15		4	5	
FIB 5X60P15		5	6	
FIB 7X60P15		7	6	

FIB.032" – 5F catheters (1,1-1,25mm) – GW .038"			
Reference	PrimarySpring Ø	Coil Ø	Coil length cm
FIB 4X50P32	0.82 mm / [.032"]	4	5
FIB 5X60P32		5	6
FIB 6X60P32		6	6
FIB 8X60P32		8	6
FIB 8X100P32		8	10
FIB 8X200P32		8	20
FIB 12X200P32		12	20
FIB 15X250P32		15	25
GW: Guide Win			



CIRRUS (with fibers)

FLOW COILS Range (Free spirales for injection)

These are soft platinum coils designed to be delivered by injection through BALT MAGIC (1,2F & 1,5F) & SONIC (1,2F & 1,5F) microcatheters.

SPIF.010" - MAGIC1,2F/1,5F or SONIC1,2F/1,5F

F	Reference	PrimarySpring Ø	Coil Ø mm	Coil length cm
S	SPIF 2,5X5P10		2,5	5
S	SPIF 2,5X10P10	0.19 mm / [.007"]	2,5	10
S	SPIF 2,5X20P10		2,5	20



10, rue de la Croix Vigneron 95160 Montmorency • France Tél. : +33 (0)1 39 89 46 41 Fax : +33 (0)1 34 17 03 46 www.balt.fr The content of this document, in particular data, information, trademarks and logos are BALT EXTRUSION SAS's (hereinafter "BALT") sole property. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT's copyrights and other intellectual proprietary rights. This document with associated pictures are non-contractual and are solely dedicated to healthcare professionals and BALT's distributors. The products commercialized by BALT shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes.