ordering information

One of the most complete coil offers on the market

	Model No.	Implant size
stige s mplex a Pack	PRES0202CXPPL	2mm x 2cm
	PRES0204CXPPL	2mm x 4cm
	PRES0206CXPPL	2mm x 6cm
	PRES0307CXPPL	3mm x 7cm
	PRES0310CXPPL	3mm x 10cm
	Model No.	Implant size
	PRES0204HXPPL	2mm x 4cm
stige	PRES0204HXPPL PRES0206HXPPL	2mm x 4cm 2mm x 6cm
s		
s ical	PRES0206HXPPL	2mm x 6cm
s ical	PRES0206HXPPL PRES0208HXPPL	2mm x 6cm 2mm x 8cm
s ical	PRES0206HXPPL PRES0208HXPPL PRES0310HXPPL	2mm x 6cm 2mm x 8cm 3mm x 10cm
s	PRES0206HXPPL PRES0208HXPPL PRES0310HXPPL	2mm x 6cm 2mm x 8cm 3mm x 10cm

Model No.	Implant size
PRES1240COM18	12mm x 40cm
PRES1447COM18	14mm x 47cm
PRES1655COM18	16mm x 55cm
PRES1860COM18	18mm x 60cm
PRES2065COM18	20mm x 65cm
PRES2465COM18	24mm x 65cm

	Model No.	Box conten
JLTRA Detachment Controller	ULTRA	1 detachment controller

Model No.	Implant size
PRES0407CPKPL	4mm x 7cm
PRES0410CPKPL	4mm x 10cm
PRES0413CPKPL	4mm x 13cm
PRES0510CPKPL	5mm x 10cm
PRES0517CPKPL	5mm x 17cm
PRES0620CPKPL	6mm x 20cm
PRES0720CPKPL	7mm x 20cm
PRES0830CPKPL	8mm x 30cm
PRES0930CPKPL	9mm x 30cm
PRES1030CPKPL	10mm x 30cm
PRES1040CPKPL	10mm x 40cm

Prestige

Complex Pack

Prestige

Plus Helical

Pack

۲

Model No.	Implant size
PRESO410HPKPL	4mm x 10cm
PRES0420HPKPL	4mm x 20cm
PRES0515HPKPL	5mm x 15cm
PRES0520HPKPL	5mm x 20cm
PRESO620HPKPL	6mm x 20cm
PRES0720HPKPL	7mm x 20cm
PRES0830HPKPL	8mm x 30cm
PRES0930HPKPL	9mm x 30cm
PRES1040HPKPL	10mm x 40cm

prestige plus volume and versatility



The Prestige Coil System is intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vessels. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/ EEC amended by 2007/47/EC). Manufactured by BALT USA LLC, 29 Parker, Irvine, CA 92618. Carefully read the instructions for use before use. First CE marking: 2021. For a list of indications, warnings, precautions and contraindications, refer to IFU-060.

indications, warnings, precoutions and contraindications, refer to IFU-060. The content of this document, in particular data, information, trademarks and logos is BALT S.A.S and affiliates' cole property. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT S.A.S and affiliates' copyrights and other intellectual proprietary rights @2021 BALT S.A.S and affiliates' all rights reserved. This document with associated partures is non-contractual and is solely dedicated to healthcare professionals and BALT S.A.S and affiliates' distributors. The products commercialized by BALT S.A.S and affiliates shall exclusively be used in accordance with the package inserts included in the boxes. DC0576B (05/2021)

Balt

10, rue de la Croix Vigneron, 95160 Montmorency, France Tél. : +33 (0)1 39 89 46 41 Fax : +33 (0)1 34 17 03 46





۲

Со

Xtro

Pre Plu He

Xtro

Complex

۲

prestige plus

Intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vessels.

Design

Premium coil softness through ariations of filar and primary wind for deliverability and conformability

Body coil/hypotube pusher design provides **increased ease of delivery**, even in tortuous anatomy

Smooth, precise embolization thanks to implant softness and enhanced pusher technology

0

Optimal detachment

Instantaneous detachment in under 1second* allows for rapid treatment & efficient case times

Reliable detachment* system provides confidence during coiling cases

Complete offer

With a full range of Complex and Helical coils to accommodate c variety of clinical application

Unique large volume coils* (Complex18 and Complex Pack) for dense packing



۲

entire size range deliverable through

.0165"-.021/.022" catheters

۲