

# Llaves de Paso de MERIT

**DEBENE**  
*La mejor calidad en alta complejidad*




**DEBENE**

# Llaves de Paso

- ❑ **Introducción**
- ❑ Llaves MANIFOLD (de 2 y 3 puertos)
- ❑ Llave de 3 Vías MARQUIS (“Stopcock”)

# Calidad Merit de Llaves



# Llaves de Paso

- ❑ Introducción
- ❑ Llaves MANIFOLD (de 2 y 3 puertos)
- ❑ Llave de 3 Vías MARQUIS (“Stopcock”)

# 502HN-R – Manifold de 2 Puertos

## 502HN-R

**50** = 500 PSI

**2** = 2 Puertos

**H** = Half Body

**N** = Manija ON

**R** = Rotación Derecha


Manija ON

Airless Rotator

Half Body  
(Medio Cuerpo)

500 psi (color  
celeste opaco)

Puerto Macho  
(otros puertos = Hembra)



# 503HN-R – Manifold de 3 Puertos

## 503HN-R

**50** = 500 PSI

**3** = 3 Puertos

**H** = Half Body

**N** = Manija ON

**R** = Rotación Derecha


Manija ON

Airless Rotator

Half Body  
(Medio Cuerpo)

500 psi (color  
celeste opaco)

Puerto Macho  
(otros puertos = Hembra)



# Llaves de Paso

- ❑ Introducción
- ❑ Llaves MANIFOLD (de 2 y 3 puertos)
- ❑ Llave de 3 Vías MARQUIS (“Stopcock”)

# MRRC – Llave de 3 Vías


## M3RRC

**M** = Presión Mediana (500 PSI)

**3** = 3 Vías

**RRC** = Right Rotating Connector

Rotador Libre de Aire  
Transparencia Conducto Interior  
Manija a forma de Ala de Tiburoón



Manija  
(Acetal)

**RRC**  
Right Rotating  
Connector  
(MALE LUER)

AIRLESS ROTATOR  
CLEAR INNER LUMEN  
SHARKS FIN HANDLE



# Otras Llaves de 3 Vías

<https://www.merit.com/wp-content/uploads/2013/11/400177018-D-Marquis-stopcock-brochure.pdf>



S1SNC — Standard pressure, one-way with swivel nut connector

M1SNC — Medium pressure, one-way with swivel nut connector

H1FLC — High pressure, one-way with fixed luer connector

H1RC — High pressure, one-way with rotating connector

S3SNC — Standard pressure, three-way with swivel nut connector

S4SNC — Standard pressure, four-way with swivel nut connector

M3SNC — Medium pressure, three-way with swivel nut connector

H3SNC — High pressure, three-way with swivel nut connector

M3RRC — Medium pressure, three-way with right rotating connector

H3RRC — High pressure three-way with right rotating connector

Other configurations are available, ask your customer service representative.

S3FLC — Standard pressure, three-way with fixed luer connector

S4FLC — Standard pressure, four-way with fixed luer connector

M3FLC — Medium pressure, three-way with fixed luer connector

H3FLC — High pressure, three-way with fixed luer connector

En stock Argentina

**DEBENE**

# Llaves de Paso

- ❑ Introducción
- ❑ Llaves MANIFOLD (de 2 y 3 puertos)
- ❑ Llave de 3 Vías MARQUIS (“Stopcock”)